

## BIOCOMPATIBILITY CERTIFICATE

Testmaterial:

Elastosil® LR3003/20 A,B Control No: 1823SL

Manufacturer:

WACKER Chemie GmbH, Germany WACKER Silicones Corporation, USA WACKER Chemicals East Asia Ltd., Japan

Studies performed:

The following studies were performed in order to determine the

biocompatiblity of the test item according to ISO 10993-1:

**CYTOTOXICITY** 

**HEMOLYSIS** 

**PYROGENICITY** 

**SENSITISATION** 

**DERMAL IRRITATION** 

**IMPLANTATION (90 days)** 

Results:

The test item showed no relevant biological effects in the studies

performed.

Suitable for use in manufacturing of medical devices for short term (no longer than 30 days) implant applications in accordance with the Wacker Health Care Letter and Guidelines.

**BSL BIOSERVICE Scientific Laboratories GmbH Munich** 

Behringstraße 6

D-82152 Planegg

Dr. Albrecht Poth

Manager Biological Safety Testing

Date: December 28, 1998